

		A		B	C
		TUB UNITS		INLET FILTERS	OUTLET FILTERS
		SIZES	MATERIAL		
TUB ASSEMBLIES	FILTER BAG OPTIONS	2 Litre Tub and Base Unit	Glass 50035006 Stainless Steel 50035005	SS 60 Mesh Support Filter (50035113) with 45 Micron Nylon Inlet Filter (50035110) as standard; Also available SS 250 Mesh 50035114 SS 500 Mesh 50035115	<p align="center">Large Filter Bags</p> Nomex 50035407 Nylon 50035400 Polypropylene 50035404 Terylene 50035402
		5 Litre Tub and Base Unit	Glass 50035009 Stainless Steel 50035008		
		Tub assembly 5 Litre with GL32 port supplied with screw fitted blank plate <i>To fit Humidity Probe 50186500</i>	Glass 50035010 (Must use Model 501 FBD to use Humidity Probe)		
		Multi Tub Unit 50035011	Stainless Steel Tubs 4 X 50035012	SS 60 Mesh Support Filter as standard	
	Glass Tubs 4 X 50035013		SS 60 Mesh Support Filter (50035309) with 45 Micron Nylon Inlet Filter (50035305) as standard; Also available SS 250 Mesh 50035307 SS 500 Mesh 50035308 3 Micron Polyester 50035306		
	SEALED TUB OPTIONS <small>(Includes Silicon O-ring & ground glass flange with clamp)</small>	5 Litre Sealed Glass Tub and Base with Clamp (50035014)		SS 60 Mesh Support Filter (50035113) with 45 Micron Nylon Inlet Filter (50035110) as standard; Also available SS 250 Mesh 50035114 SS 500 Mesh 50035115	Top Cap with 45 Micron Nylon Filter 50035020 Top Cap with SS 250 Mesh Filter 50035021 Top Cap with SS 500 Mesh Filter 50035022 Top Cap with Nylon Filter Bag 50035023 Top Cap + 3 Micron Polyester Filter 50035024
		Tub Assembly 5 Litre Glass Moisture Analysis (Includes ports for sample and Humidity Probe. Both ports supplied with screw fitted blank plates)	Glass 50135020 (Must use Model 501 FBD to use Humidity Probe <u>50186500</u>)	3 Micron Polyester Filter for smallest particle size 50035120	
		Low Density Classifier Assembly	Glass Tub 50035049		
		Multi Tub Unit 50035011	Glass Tub 4 x 50035033	Mini Tub, Sealed Tub, Inlet and Outlet sets Mini 45 Micron Nylon Filter Set 50035310 Mini 3 Micron Polyester Filter Set 50035311 Mini SS 250 Mesh Filter Set 50035312 Mini SS 500 Mesh Filter Set 50035313 (SS 60 Mesh Support Filter (50035309) and 45 Micron Nylon Inlet Filter (50035305) as standard)	